

OPERATION TABLES FOR ADVANCED SURGERIES

MOTORIZED OPERATION TABLE - ELECTRO WITH MANUAL OVER RIDE

OPTAB EMO 01



- 3 way Operated Operation Table :
 - ✓ Electrically
 - ✓ On Battery Back up
 - ✓ Manually (With Help of Crank)
- Change of operating mode by selector switch on failure of electric or functioning of remote.

Anaesthetic Screen	1 PC
Arm Rest	1 Pair
Shoulder Support with pad	1 Pair
Lateral Support with pad	1 Pair
Waterproof Mattress	1 Set
Lithotomy Crutches	1 Pair
Full body covered by S S Sheet	
Adjustable, Removable Interchangeable Head & Foot Top	
With castor Wheel for easy mobility along with stable floor locking	

MOTORIZED OPERATION TABLE - ELECTRO

OPTAB ELECTRO 01



Salient Features :

- Fully Electric with Cable Remote
- Off line Battery Back up for emergency
- Weight Capacity up to 200 Kg
- Four Separate Actuators for
 - ✓ Vertical movement
 - ✓ Trendelburg & Reverse Trend.
 - ✓ Side Tilting
 - ✓ Chair Position
- X-ray permeable Table Top without any cross bars.
- Eccentric Pillar to make the table real IITV Compatible
- Most convenient for all type of major & minor surgery

Technical Parameters

L : 1980mm	W : 530mm
H : 760 - 1010mm / 660 - 910mm	
Lateral Tilt : $\pm 20^\circ$	
Trend. / R Trend. : $\pm 25^\circ$	
Back Section: $+85^\circ / -40^\circ$	
Head & Leg : $\pm 90^\circ$	

On Demand Available :

- ✓ ACP Control or additional Remote
- ✓ Zero Position facility
- ✓ Top Slide Facility

ACCUMAX SURGICALS

Manufacturer of : HOSPITAL EQUIPMENTS & FURNITURE



Admin & Works :

55, Amar Estate, Opp. Forge & Blower Comp.
Nr. Memco, Naroda Road, Ahmedabad-25. INDIA.
Ph : +91-79-22203485

Works & Display :

C/127, Vishala Ind. Estate,
Outside Kathawada G.I.D.C., S.P Ring Road,
Odhav, Ahmedabad-342415. Ph : +91-7929092177

Mo.: +91 99798 58530 | E-mail : accumaxsurgicals@hotmail.com | Website : www.accumaxsurgicals.com