

STRETCHER ON TROLLEY

New Design – Folding



- Compaq with Excellent Look
- Made of rectangular pipe for better strength
- **Non detachable** S.S. Top with folding flaps
- 6" wheel castors

Regular



- M S Powder Coated with S.S. Handle
- Detachable Top
- Overall size : 2030mm x 600 mm x 810 mm
- Wheel Castors : 4" / 6" / 8" as per confirmation
- Optionally Available : Superimposed S.S. Top

Regular – Folding



- M S Powder Coated with S S Handle
- Detachable Top
- Top with Folding sides to minimize storage size.
- Overall size : 2030mm x 600 mm x 810 mm
- Wheel Castors : 4" / 6" / 8" as per confirmation
- Optionally Available : Superimposed S S Top

ACCUMAX SURGICALS

Manufacturer & Supplier of : HOSPITAL EQUIPMENTS & FURNITURE



AN ISO 9001:2008 & ISO 13485:2012 Certified Company



ACCUMAX

Admin & Works :

55, Amar Estate, Opp. Forge & Blower Comp,
Nr. Memco, Naroda Road, Ahmedabad-25. INDIA.
Ph : +91-79-22203485 M. : 9979858530

Works & Display :

C/127, Vishala Industrial Estate(E),
Outside Kathawada G.I.D.C., S.P. Ring Road,
Odhav, Ahmedabad-342415. Ph : +91-79-29092177

E-mail : accumaxsurgicals@hotmail.com

Website : www.accumaxsurgicals.com

Products and specifications are Subject to change without notice