

STRETCHER ON TROLLEY

New Design – Folding - STR ND FLD 01



- Compaq with Excellent Look
- Made of rectangular pipe for better strength
- **Non detachable** S.S. Top with folding flaps
- 6" wheel castors

Regular - STR REG 01



- M S Powder Coated with S.S. Handle
- Detachable Top
- Overall size : 2030mm x 600 mm x 810 mm
- Wheel Castors : 4" / 6" / 8" as per confirmation
- Optionally Available : Superimposed S.S. Top

Regular – Folding - STR REG FLD 01



- M S Powder Coated with S S Handle
- Detachable Top
- Top with Folding sides to minimize storage size.
- Overall size : 2030mm x 600 mm x 810 mm
- Wheel Castors : 4" / 6" / 8" as per confirmation
- Optionally Available : Superimposed S S Top

ACCUMAX SURGICALS

Manufacturer of : HOSPITAL EQUIPMENTS & FURNITURE



ACCUMAX

Admin & Works :

55, Amar Estate, Opp. Forge & Blower Comp.
Nr. Memco, Naroda Road, Ahmedabad-25. INDIA.
Ph : +91-79-22203485

Works & Display :

C/127, Vishala Ind. Estate,
Outside Kathawada G.I.D.C., S.P Ring Road,
Odhav, Ahmedabad-342415. Ph : +91-7929092177

Mo.: +91 99798 58530 | E-mail : accumaxsurgicals@hotmail.com | Website : www.accumaxsurgicals.com