

# EXAMINATION TABLE

## GYNEC EXAMINATION TABLE

GY EXT 01



- M.S. Powder Coated Plain Structure
- Over all Size : 72" x 20" x 32"
- Head Portion of 18" length adjustable on Ratchet
- Optionally Available : Superimposed S.S. Top & Mattress

## PLAIN EXAMINATION TABLE

PL EXT 01



- M.S. Powder Coated Structure
- Perineal cut on leg Side
- Over all Size : 72" x 24" x 32"
- Head Portion of 18" length adjustable on Ratchet
- Optionally Available : Superimposed S.S. Top & Mattress

## ACCUMAX SURGICALS

Manufacturer of : HOSPITAL EQUIPMENTS & FURNITURE



### Admin & Works :

55, Amar Estate, Opp. Forge & Blower Comp,  
Nr. Memco, Naroda Road, Ahmedabad-25. INDIA.  
Ph : +91-79-22203485

### Works & Display :

C/127, Vishala Ind. Estate,  
Outside Kathawada G.I.D.C., S.P Ring Road,  
Odhav, Ahmedabad-342415. Ph : +91-7929092177

Mo.: +91 99798 58530 | E-mail : [accumaxsurgical@hotmail.com](mailto:accumaxsurgical@hotmail.com) | Website : [www.accumaxsurgical.com](http://www.accumaxsurgical.com)