

ICU, RECOVERY & FOWLER BED

ICU BED (Without Hi-Lo) ICU BD 01



- M S Powder Coated Structure with Four Section Perforated Top
- Three Separate Screws for Back Rest , Knee Rest & Trendlenburg - Reverse Trendlenburg
- Four Provisions for Saline Stand on Main Frame
- Over all size : 84" x 36" x 24" (H) **without Hi – Lo Facility**
- 5" Swiveling Castors : Regular / P U (Optional)
- Available with **M S Bow, S.S. Bow – Lam Panel, Fully S.S. Bow or ABS Panel** as per Order
- Optional Supply : Side Railing, On Bed Saline Stand , Mattress

RECOVERY BED RCV BD 01



- M S Powder Coated Structure with Two Section Perforated Top
- Two Separate Screws for Back Rest Trendlenburg - Reverse Trendlenburg
- Four Provisions for Saline Stand on Main Frame
- Over all size : 78 " x 36" x 24" (H)
- 5" Swiveling Castors : Regular / P U (Optional)
- Available with **M S Bow, S.S. Bow – Lam Panel, Fully S.S. Bow or ABS Panel** as per Order
- Optional Supply : Side Railing, On Bed Saline Stand, Mattress

FOWLER BED FWL BD 01



- M S Powder Coated Structure with Four Section Perforated Top
- Two Separate Screws for Back Rest & Knee Rest
- Four Provisions for Saline Stand on Main Frame
- Over all size : 84" x 36" x 24" (H)
- Available with **M S Bow, S.S. Bow – Lam Panel, Fully S.S. Bow or ABS Panel** as per Order
- Optional Supply : 5" Swiveling Castors (Regular or PU), Side Railing, On Bed Saline Stand , Mattress



ABS PNL 01

ABS Panel



SS BLM PNL 01

S.S. Bow-Lam Panel



FUL SS BW 01

Fully S.S. Bow



MS BW 01

M.S. Bow

ACCUMAX SURGICALS

Manufacturer of : HOSPITAL EQUIPMENTS & FURNITURE



ACCUMAX

Admin & Works :

55, Amar Estate, Opp. Forge & Blower Comp,
Nr. Memco, Naroda Road, Ahmedabad-25. INDIA.
Ph : +91-79-22203485

Works & Display :

C/127, Vishala Ind. Estate,
Outside Kathawada G.I.D.C., S.P Ring Road,
Odhav, Ahmedabad-342415. Ph : +91-7929092177

Mo.: +91 99798 58530 | E-mail : accumaxsurgicals@hotmail.com | Website : www.accumaxsurgicals.com