

ICU, RECOVERY & FOWLER BED

ICU Bed (Without Hi-Lo)



- M S Powder Coated Structure with Four Section Perforated Top
- Three Separate Screws for Back Rest , Knee Rest & Trendlenburg - Reverse Trendelnburg
- Four Provisions for Saline Stand on Main Frame
- Over all size : 84" x 36" x 24" (H) **without Hi – Lo Facility**
- 5" Swiveling Castors : Regular / P U (Optional)
- Available with **M S Bow, S.S. Bow – Lam Panel, Fully S.S. Bow or ABS Panel** as per Order
- Optional Supply : Side Railing, On Bed Saline Stand , Mattress

Recovery Bed



- M S Powder Coated Structure with Two Section Perforated Top
- Two Separate Screws for Back Rest Trendlenburg - Reverse Trendelnburg
- Four Provisions for Saline Stand on Main Frame
- Over all size : 78 " x 36" x 24" (H)
- 5" Swiveling Castors : Regular / P U (Optional)
- Available with **M S Bow, S.S. Bow – Lam Panel, Fully S.S. Bow or ABS Panel** as per Order
- Optional Supply : Side Railing, On Bed Saline Stand, Mattress

Fowler Bed



- M S Powder Coated Structure with Four Section Perforated Top
- Two Separate Screws for Back Rest & Knee Rest
- Four Provisions for Saline Stand on Main Frame
- Over all size : 84" x 36" x 24" (H)
- Available with **M S Bow, S.S. Bow – Lam Panel, Fully S.S. Bow or ABS Panel** as per Order
- Optional Supply : 5" Swiveling Castors (Regular or PU), Side Railing, On Bed Saline Stand , Mattress



ABS Panel



S.S. Bow-Lam Panel



Fully S.S. Bow



M.S. Bow

ACCUMAX SURGICALS

Manufacturer & Supplier of : HOSPITAL EQUIPMENTS & FURNITURE



AN ISO 9001:2008 & ISO 13485:2012 Certified Company

Admin & Works :

55, Amar Estate, Opp. Forge & Blower Comp,
Nr. Memco, Naroda Road, Ahmedabad-25. INDIA.
Ph : +91-79-22203485 M. : 9979858530

Works & Display :

C/127, Vishala Industrial Estate(E),
Outside Kathawada G.I.D.C., S.P. Ring Road,
Odhav, Ahmedabad-342415. Ph : +91-79-29092177



ACCUMAX

E-mail : accumaxsurgicals@hotmail.com

Website : www.accumaxsurgicals.com

Products and specifications are Subject to change without notice