

# SEMIFOWLER, PLAIN & PEDIATRIC BED

## SEMIFOWLER BED SMFWL BD 01



- M S Powder Coated Structure with Two Section Perforated Top
- Adjustment of Back Rest up to 85 Deg by means of Screw
- Four Provisions for Saline Stand on Main Frame
- Over all size : 78" x 36" x 24" (H)
- Available with **M S Bow, S.S. Bow – Lam Panel, Fully S.S. Bow or ABS Panel** as per Order
- Optional Supply : 5" Swiveling Castors (Regular or PU), Side Railing, On Bed Saline Stand, Mattress

## PLAIN SHEET METAL BED PLSM BD 01



- M S Powder Coated Structure with Plain Perforated Top
- Four Provisions for Saline Stand on Main Frame
- Over all size : 78" x 36" x 24" (H)
- Available with **M S Bow, S.S. Bow – Lam Panel, Fully S.S. Bow or ABS Panel** as per Order
- Optional Supply : 5" Swiveling Castors (Regular or PU), Side Railing, On Bed Saline Stand, Mattress



**ABS PNL 01**

**ABS Panel**



**SS BLM PNL 01**

**S.S. Bow-Lam Panel**



**FUL SS BW 01**

**Fully S.S. Bow**



**MS BW 01**

**M.S. Bow**

## PEDIATRIC BED WITH RAILING PED BD 01



- M S Powder Coated Structure with Plain Perforated Top
- Drop down railing on Two Sides
- Four Provisions for Saline Stand on Main Frame
- Over all size : 63" x 30" x 24" (H)
- Optional Supply : On Bed Saline Stand, Mattress

## ACCUMAX SURGICALS

Manufacturer of : HOSPITAL EQUIPMENTS & FURNITURE



**ACCUMAX**

### Admin & Works :

55, Amar Estate, Opp. Forge & Blower Comp,  
Nr. Memco, Naroda Road, Ahmedabad-25. INDIA.  
Ph : +91-79-22203485

### Works & Display :

C/127, Vishala Ind. Estate,  
Outside Kathawada G.I.D.C., S.P Ring Road,  
Odhav, Ahmedabad-342415. Ph : +91-7929092177

Mo.: +91 99798 58530 | E-mail : [accumaxsurgicals@hotmail.com](mailto:accumaxsurgicals@hotmail.com) | Website : [www.accumaxsurgicals.com](http://www.accumaxsurgicals.com)