

BED SIDE UTILITIES

Bed Side Locker



Deluxe

M S Powder Coated Structure with S S Top
One Drawer on upper & Cabinet on lower portion
Mounted on 50 mm Swiveling Castors
Over all size : 15" x 15" x 32"



Regular

M S Powder Coated Structure with S S Top
Cabinet on middle portion
Over all size : 17" x 15" x 32"



Open – Fully S.S.

Fully S S Structure with Three Shelves
Special for Open & visible cleanliness
Over all size : 15" x 12" x 32"

Attendant Bed



- M S Powder Coated Structure with 2" Mattress
- Comfortably useful for visitor as well as attendant
- Over all Size : 72" x 24" x 18"

Over Bed Table-Gear Model



- Height adjustment by means of lever through M S Bevel- Gear Set
- Framework of M S P C Structure mounted on 2" wheel
- Size : 32" x 14" Laminated Top with 24" x 14" with angular adjustment

Bedside Stool



- M S Powder Coated
- Over all size : 12" x 12" x 20" (H)
- Optionally available Superimposed S S Top

ACCUMAX SURGICALS

Manufacturer & Supplier of : HOSPITAL EQUIPMENTS & FURNITURE



AN ISO 9001:2008 & ISO 13485:2012 Certified Company



ACCUMAX

Admin & Works :

55, Amar Estate, Opp. Forge & Blower Comp,
Nr. Memco, Naroda Road, Ahmedabad-25. INDIA.
Ph : +91-79-22203485 M. : 9979858530

Works & Display :

C/127, Vishala Industrial Estate(E),
Outside Kathawada G.I.D.C., S.P. Ring Road,
Odhav, Ahmedabad-342415. Ph : +91-79-29092177

E-mail : accumaxsurgicals@hotmail.com

Website : www.accumaxsurgicals.com

Products and specifications are Subject to change without notice