

BED SIDE UTILITIES

BED SIDE LOCKER



Deluxe - DEL LCR 01

(Size : 15" x 15" x 32")

M S Powder Coated Structure with S S Top
One Drawer on upper & Cabinet on lower portion
Mounted on 50 mm Swiveling Castors
Over all size : 15" x 15" x 32"



Regular - REG LCR 01

(Size : 17" x 15" x 32")

M S Powder Coated Structure with S S Top
Cabinet on middle portion
Over all size : 17" x 15" x 32"



Open-Fully S.S. - OP SS LCR 01

(Size : 15" x 12" x 32")

Fully S S Structure with Three Shelves
Special for Open & visible cleanliness
Over all size : 15" x 12" x 32"

ATTENDANT BED - ATD BD 01



- M S Powder Coated Structure with 2" Mattress
- Comfortably useful for visitor as well as attendant
- Over all Size : 72" x 24" x 18"

OVER BED TABLE-GEAR MODEL

OVBD WOOD 01

OVBD ABS 01



- Height adjustment by means of lever through M S Bevel- Gear Set
- Framework of M S P C Structure mounted on 2" wheel
- Size : 32" x 14" Laminated Top with 24" x 14" with angular adjustment

BEDSIDE STOOL BD SD STL 01



- M S Powder Coated
- Over all size : 12" x 12" x 20" (H)
- Optionally available Superimposed S S Top

ACCUMAX SURGICALS

Manufacturer of : HOSPITAL EQUIPMENTS & FURNITURE

Admin & Works :

55, Amar Estate, Opp. Forge & Blower Comp,
Nr. Memco, Naroda Road, Ahmedabad-25. INDIA.
Ph : +91-79-22203485

Works & Display :

C/127, Vishala Ind. Estate,
Outside Kathawada G.I.D.C., S.P Ring Road,
Odhav, Ahmedabad-342415. Ph : +91-7929092177



Mo.: +91 99798 58530 | E-mail : accumaxsurgicalsh@hotmail.com | Website : www.accumaxsurgicalsh.com