

BED SIDE SCREEN & SALINE STAND

BED SIDE SCREEN - BD SD SCR N 01



- Steel tubular structure made of 19 mm diameter steel tubes
- Mounted on 2" Wheel Castors
- Three powder coated sections provided with spring wires & hooks
- Over all size : 2440 mm x 1675 mm

SALINE STAND



Fully S.S. SLN SS 01

- Fully S S Structure with S S Rod
- Preferable for O T & Deluxe, Sp. or Semi Sp. Rooms
- 5 Nos. of 2" Swiveling Castors Or Buffer

Fully S.S.

S.S. with
Fiber Base

M S P C
with Castors

Regular

S.S. with Fiber Base SLN SS FB 01

- Stainless Steel Pipe & Rod with Fiber base
- Ensure safety for heavy bottle due to wide base
- 5 Nos. of 2" Swiveling Castors



MSPC with Castors SLN SS MSPC 01

- M S Powder Coated Structure with S S Rod
- Stronger with less Weight
- 5 Nos. of 2" Swiveling Castors

Regular SLN REG 01

- M S Powder Coated Pipe with M S Chrome Plated Rod
- Cast Iron Round Base



ACCUMAX SURGICALS

Manufacturer of : HOSPITAL EQUIPMENTS & FURNITURE



Admin & Works :
55, Amar Estate, Opp. Forge & Blower Comp,
Nr. Memco, Naroda Road, Ahmedabad-25. INDIA.
Ph : +91-79-22203485

Works & Display :
C/127, Vishala Ind. Estate,
Outside Kathawada G.I.D.C., S.P Ring Road,
Odhav, Ahmedabad-342415. Ph : +91-7929092177

Mo.: +91 99798 58530 | E-mail : accumaxsurgicals@hotmail.com | Website : www.accumaxsurgicals.com