

# BED SIDE SCREEN & SALINE STAND

## Bed Side Screen



- Steel tubular structure made of 19 mm diameter steel tubes
- Mounted on 2" Wheel Castors
- Three powder coated sections provided with spring wires & hooks
- Over all size : 2440 mm x 1675 mm

## Saline Stand

### Fully S.S.

- Fully S S Structure with S S Rod
- Preferable for O T & Deluxe, Sp. or Semi Sp. Rooms
- 5 Nos. of 2" Swiveling Castors Or Buffer



### S.S. with Fiber Base

- Stainless Steel Pipe & Rod with Fiber base
- Ensure safety for heavy bottle due to wide base
- 5 Nos. of 2" Swiveling Castors



### MSPC with Castors

- M S Powder Coated Structure with S S Rod
- Stronger with less Weight
- 5 Nos. of 2" Swiveling Castors



### Regular

- M S Powder Coated Pipe with M S Chrome Plated Rod
- Cast Iron Round Base



## ACCUMAX SURGICALS

Manufacturer & Supplier of : HOSPITAL EQUIPMENTS & FURNITURE



AN ISO 9001:2008 & ISO 13485:2012 Certified Company



## ACCUMAX

### Admin & Works :

55, Amar Estate, Opp. Forge & Blower Comp,  
Nr. Memco, Naroda Road, Ahmedabad-25. INDIA.  
Ph : +91-79-22203485 M. : 9979858530

### Works & Display :

C/127, Vishala Industrial Estate(E),  
Outside Kathawada G.I.D.C., S.P. Ring Road,  
Odhav, Ahmedabad-342415. Ph : +91-79-29092177

E-mail : [accumaxsurgicals@hotmail.com](mailto:accumaxsurgicals@hotmail.com)

Website : [www.accumaxsurgicals.com](http://www.accumaxsurgicals.com)

Products and specifications are Subject to change without notice