

HYDRAULIC OPERATION TABLE

HYDRAULIC OPERATION TABLE : "C" ARM - OPTAB CARM 01



Technical Parameters

L : 2100mm	W : 500mm
H : 750 - 1050mm / 710 - 960mm	
Lateral Tilt : $\pm 20^\circ$	
Trend. / R Trend. : $\pm 25^\circ$	
Back Section: $+85^\circ / -40^\circ$	
Head & Leg : $\pm 90^\circ$	

Salient Features :

- Manually Operated smoothly & jerk free.
- Eccentric Pillar to make the table real IITV Compatible.
- X-ray permeable Table Top without any cross bars
- Most convenient for all type of general & uro surgery
- Weight Capacity up to 200 Kg
- Manually operated :
 - ✓ Vertical movement by foot padde
 - ✓ Trendelburg & Reverse Trendelburg
 - ✓ Side Tilting
 - ✓ Chair Position

HYDRAULIC OPERATION TABLE : SUPER DELUXE

OPTAB SUPDL 01



Technical Parameters

L : 1830mm	W : 510mm
H : 750 - 1050mm / 660 - 910mm	
Lateral Tilt : $\pm 20^\circ$	
Trend. / R Trend. : $\pm 25^\circ$	
Back Section: $+85^\circ / -40^\circ$	
Head & Leg : $\pm 90^\circ$	

Salient Features :

- Fully manual with Hydraulic Pump and Smooth uniform gear mechanism
- Best suited Economy table for all non uro surgery
- Elevated by single lever double action horizontal hydraulic pump
- Smooth individual side controls for :
 - ✓ Trendelburg /Reverse trendelburg
 - ✓ Flex/ reflex
 - ✓ Chair position
 - ✓ Side Tilting

ACCUMAX SURGICALS

Manufacturer of : HOSPITAL EQUIPMENTS & FURNITURE



Admin & Works :

55, Amar Estate, Opp. Forge & Blower Comp,
Nr. Memco, Naroda Road, Ahmedabad-25. INDIA.
Ph : +91-79-22203485

Works & Display :

C/127, Vishala Ind. Estate,
Outside Kathawada G.I.D.C., S.P Ring Road,
Odhav, Ahmedabad-342415. Ph : +91-7929092177

Mo.: +91 99798 58530 | E-mail : accumaxsurgicals@hotmail.com | Website : www.accumaxsurgicals.com