

ORTHOPAEDIC OPERATION TABLE

Hydraulic : "C" Arm Table with Orthopedic Attachment




Salient Features :

- Manually Operated uniformly & jerk free.
- Adjustable, Removable as well as interchangeable Head & Foot Top as per requirement
- Single lever floor locking
- Full body covered by S.S. Sheet
- Weight Capacity up to 200 Kg
- Vertical movement by Hydraulic Pump and other by Crank.
- Eccentric Pillar to make the table Compatible with Orthopedic Attachment.
- Fully S.S. Orthopedic Attachment compatible with IITV.
- logically designed Orthopaedic Attachment to suit all procedures

Technical Parameters

L : 2100mm	W : 500mm
H : 750 - 1050mm / 710 - 960mm	
Lateral Tilt : $\pm 20^\circ$	
Trend. / R Trend. : $\pm 25^\circ$	
Back Section: $+85^\circ / -40^\circ$	
Head & Leg : $\pm 90^\circ$	



Standard Orthopaedic Accessories

- Radiolucent Pelvic Post with centric & eccentric positioning
- Standard Leg Traction Unit – Telescopic type (for Tibia & Femur Traction)
- Popliteal Support for Tibia Close nailing
- S T Pin Holding clamp
- X-ray permeable Large size Arm Rest.
- Attachment for Elbow Surgery



ACCUMAX SURGICALS

Manufacturer & Supplier of : HOSPITAL EQUIPMENTS & FURNITURE



AN ISO 9001:2008 & ISO 13485:2012 Certified Company



ACCUMAX

Admin & Works :

55, Amar Estate, Opp. Forge & Blower Comp,
Nr. Memco, Naroda Road, Ahmedabad-25. INDIA.
Ph : +91-79-22203485 M. : 9979858530

Works & Display :

C/127, Vishala Industrial Estate(E),
Outside Kathawada G.I.D.C., S.P. Ring Road,
Odhav, Ahmedabad-342415. Ph : +91-79-29092177

E-mail : accumaxsurgicals@hotmail.com

Website : www.accumaxsurgicals.com

Products and specifications are Subject to change without notice