

# ORTHOPAEDIC OPERATION TABLE

## ORTHOPAEDIC OPERATION TABLE - OPTAB FRT 01




### Salient Features :

- Manually Operated uniformly & jerk free.
- Adjustable, Removable as well as interchangeable Head & Foot Top as per requirement
- Single lever floor locking
- Full body covered by S.S. Sheet
- Weight Capacity up to 200 Kg
- Vertical movement by Hydraulic Pump and other by Crank.
- Eccentric Pillar to make the table Compatible with Orthopedic Attachment.
- Fully S.S. Orthopedic Attachment compatible with IITV.
- logically designed Orthopaedic Attachment to suit all procedures

### Technical Parameters

L : 2100mm	W : 500mm
H : 750 - 1050mm / 710 - 960mm	
Lateral Tilt : $\pm 20^\circ$	
Trend. / R Trend. : $\pm 25^\circ$	
Back Section: $+85^\circ / -40^\circ$	
Head & Leg : $\pm 90^\circ$	



### Standard Orthopaedic Accessories

- Radiolucent Pelvic Post with centric & eccentric positioning
- Standard Leg Traction Unit – Telescopic type (for Tibia & Femur Traction)
- Popliteal Support for Tibia Close nailing
- S T Pin Holding clamp
- X-ray permeable Large size Arm Rest.
- Attachment for Elbow Surgery



## ACCUMAX SURGICALS

Manufacturer of : HOSPITAL EQUIPMENTS & FURNITURE



### Admin & Works :

55, Amar Estate, Opp. Forge & Blower Comp,  
Nr. Memco, Naroda Road, Ahmedabad-25. INDIA.  
Ph : +91-79-22203485

### Works & Display :

C/127, Vishala Ind. Estate,  
Outside Kathawada G.I.D.C., S.P Ring Road,  
Odhav, Ahmedabad-342415. Ph : +91-7929092177

Mo.: +91 99798 58530 | E-mail : [accumaxsurgicals@hotmail.com](mailto:accumaxsurgicals@hotmail.com) | Website : [www.accumaxsurgicals.com](http://www.accumaxsurgicals.com)