

SUCTION MACHINE & AUTOCLAVE STERILIZER

Suction Machine



Mono

- Best Quality ¼ HP Single Phase Electric Motor
- Twin Bottle (2/1.5 Ltr) with overflow valve
- Self sealing rubber bush with nontoxic Rubber pipe
- Oil immersed Vacuum Pump
- Heavy duty castors for easy mobility



Highvac

Autoclave



12" x 12"
(Single Drum)

12" x 21"
(Double Drum)

12" x 27"
(2 ½ Drum)

12" x 32"
(Three Drum)

- Single S.S. Chamber with Gun Metal Lid
- Brass Wing nut with S.S. Stud for the best tightening
- Dual Scale Pressure Gauge.
- Weight calibrated Steam release valve for friendly working
- Additional Safety & Release Valves for safe Operation
- Available : Electric / Non Electric Model

ACCUMAX SURGICALS

Manufacturer & Supplier of : HOSPITAL EQUIPMENTS & FURNITURE



AN ISO 9001:2008 & ISO 13485:2012 Certified Company



ACCUMAX

Admin & Works :

55, Amar Estate, Opp. Forge & Blower Comp,
Nr. Memco, Naroda Road, Ahmedabad-25. INDIA.
Ph : +91-79-22203485 M. : 9979858530

Works & Display :

C/127, Vishala Industrial Estate(E),
Outside Kathawada G.I.D.C., S.P. Ring Road,
Odhav, Ahmedabad-342415. Ph : +91-79-29092177

E-mail : accumaxsurgicals@hotmail.com

Website : www.accumaxsurgicals.com

Products and specifications are Subject to change without notice