

SUCTION MACHINE & AUTOCLAVE STERILIZER

SUCTION MACHINE



- Best Quality ¼ HP Single Phase Electric Motor
- Twin Bottle (2/1.5 Ltr) with overflow valve
- Self sealing rubber bush with nontoxic Rubber pipe
- Oil immersed Vacuum Pump
- Heavy duty castors for easy mobility



HIGHVAC - SUC HVAC 01

HIGHVAC - SUC HVAC 01

AUTOCLAVE



12" x 12"
(Single Drum)
ATCL 12

12" x 21"
(Double Drum)
ATCL 21

12" x 27"
(2 ½ Drum)
ATCL 27

12" x 32"
(Three Drum)
ATCL 32

- Single S.S. Chamber with Gun Metal Lid
- Brass Wing nut with S.S. Stud for the best tightening
- Dual Scale Pressure Gauge.
- Weight calibrated Steam release valve for friendly working
- Additional Safety & Release Valves for safe Operation
- Available : Electric / Non Electric Model

ACCUMAX SURGICALS

Manufacturer of : HOSPITAL EQUIPMENTS & FURNITURE



Admin & Works :
55, Amar Estate, Opp. Forge & Blower Comp,
Nr. Memco, Naroda Road, Ahmedabad-25. INDIA.
Ph : +91-79-22203485

Works & Display :
C/127, Vishala Ind. Estate,
Outside Kathawada G.I.D.C., S.P Ring Road,
Odhav, Ahmedabad-342415. Ph : +91-7929092177

Mo.: +91 99798 58530 | E-mail : accumaxsurgicals@hotmail.com | Website : www.accumaxsurgicals.com