

# MONITOR TROLLEY

## Laparoscopy Trolley (Monitor Trolley)

Exclusive



- Logically designed with visual Height and comfortable size
- Monitor Tray Height : 54 – 56”
- Shelf Size : 24” x 18”
- Powder Coated with two excellent colors
- Four shelves with one lockable drawer – additional Shelf
- Net Floor coverage : 28” x 28”
- Electric Panel with 6 Nos. of 15/5 Amp. ISI marked Plugs–Switches on the side of Unit
- 5” Swiveling Castors, two of them with brake
- Provision for holding CQ Cylinder on the back of Unit ensuring Stable mobility.

Economy



## ECG-Multipara Trolley



- Z shaped Trolley with Excellent Look
- Three shelves with One drawer below top shelf
- 2” Swiveling castors with easy mobility
- Over all Size : 20” x 15” x 32”

**ACCUMAX SURGICALS**

Manufacturer & Supplier of : HOSPITAL EQUIPMENTS & FURNITURE



AN ISO 9001:2008 & ISO 13485:2012 Certified Company



**ACCUMAX**

**Admin & Works :**

55, Amar Estate, Opp. Forge & Blower Comp,  
Nr. Memco, Naroda Road, Ahmedabad-25. INDIA.  
Ph : +91-79-22203485 M. : 9979858530

**Works & Display :**

C/127, Vishala Industrial Estate(E),  
Outside Kathawada G.I.D.C., S.P. Ring Road,  
Odhav, Ahmedabad-342415. Ph : +91-79-29092177

E-mail : [accumaxsurgicals@hotmail.com](mailto:accumaxsurgicals@hotmail.com)

Website : [www.accumaxsurgicals.com](http://www.accumaxsurgicals.com)

Products and specifications are Subject to change without notice