

# CRASH CART & DRUG TROLLEY

## Troma Care Crash Cart (Emergency Equipment Cart)



M.S. Structure



S.S. Structure

- S.S. / M.S. tubular Framework – as per requirement
- Six Nos. of removable multicolor bins
- Two lockable separate plastic storage units each having 3 drawers
- Top drawers having containers of different size.
- 5" / 4" (2+2) Swiveling castors
- Pull out cardiac massage board above drawers
- Having S S saline Rod and M S P C Oxygen Cylinder Holder on one side
- Over all size : 950 mm x 500 mm x 1540 mm +/- 5%

## Medicine-Drug Injection Trolley



- Two lockable separate plastic storage units each having 3 drawers
- M.S. Tubular framework with two S.S. trays (Available in S.S. on demand)
- 4" Swiveling castors.
- Over all Size :

**ACCUMAX SURGICALS**

Manufacturer & Supplier of : HOSPITAL EQUIPMENTS & FURNITURE



AN ISO 9001:2008 & ISO 13485:2012 Certified Company



**ACCUMAX**

### Admin & Works :

55, Amar Estate, Opp. Forge & Blower Comp,  
Nr. Memco, Naroda Road, Ahmedabad-25. INDIA.  
Ph : +91-79-22203485 M. : 9979858530

### Works & Display :

C/127, Vishala Industrial Estate(E),  
Outside Kathawada G.I.D.C., S.P. Ring Road,  
Odhav, Ahmedabad-342415. Ph : +91-79-29092177

E-mail : [accumaxsurgicals@hotmail.com](mailto:accumaxsurgicals@hotmail.com)

Website : [www.accumaxsurgicals.com](http://www.accumaxsurgicals.com)

Products and specifications are Subject to change without notice