

# CRASH CART & DRUG TROLLEY

## TROMA CARE CRASH CART (EMERGENCY EQUIPMENT CART)



S.S. Structure  
CRCT SS 01



M.S. Structure  
CRCT MS 01

- S.S. / M.S. tubular Framework – as per requirement
- Six Nos. of removable multicolor bins
- Two lockable separate plastic storage units each having 3 drawers
- Top drawers having containers of different size.
- 5" / 4" (2+2) Swiveling castors
- Pull out cardiac massage board above drawers
- Having S S saline Rod and M S P C Oxygen Cylinder Holder on one side
- Over all size : 950 mm x 500 mm x 1540 mm +/- 5%

## MEDICINE-DRUG INJECTION TROLLEY - MDI TR MS 01



- Two lockable separate plastic storage units each having 3 drawers
- M.S. Tubular framework with two S.S. trays (Available in S.S. on demand)
- 4" Swiveling castors.
- Over all Size :

## ACCUMAX SURGICALS

Manufacturer of : HOSPITAL EQUIPMENTS & FURNITURE



**Admin & Works :**  
55, Amar Estate, Opp. Forge & Blower Comp,  
Nr. Memco, Naroda Road, Ahmedabad-25. INDIA.  
Ph : +91-79-22203485

**Works & Display :**  
C/127, Vishala Ind. Estate,  
Outside Kathawada G.I.D.C., S.P Ring Road,  
Odhav, Ahmedabad-342415. Ph : +91-7929092177

Mo.: +91 99798 58530 | E-mail : [accumaxsurgicals@hotmail.com](mailto:accumaxsurgicals@hotmail.com) | Website : [www.accumaxsurgicals.com](http://www.accumaxsurgicals.com)