

LABOR TABLE

Two Section Labor Table - Fully SS / MS Powder Coated with Superimposed SS Top



- Available in Fully S.S. Structure as well as M S P C with Superimposed S.S. Top
- New improved design with Superb finish and excellent workmanship
- Whole table top capable of being tilted by suitable screw mechanism
- Foot section capable of being folded underneath to suit Lithotomy position
- Back Section adjustable on rack
- Protective side rails on three sides
- Full Size : 72" x 30" x 30"
- Available with Mattress and Padded Leg Rest (Only in Fully S.S., Not in S.S. Top Model)

Plain Labor Table



- M.S. Powder Coated Structure
- M S C P Lithotomy Rod with Straps with Perineal cut on leg Side
- Over all Size : 72" x 30" x 32"
- Optionally Available : Superimposed S.S. Top & Mattress

ACCUMAX SURGICALS

Manufacturer & Supplier of : HOSPITAL EQUIPMENTS & FURNITURE



AN ISO 9001:2008 & ISO 13485:2012 Certified Company



ACCUMAX

Admin & Works :

55, Amar Estate, Opp. Forge & Blower Comp,
Nr. Memco, Naroda Road, Ahmedabad-25. INDIA.
Ph : +91-79-22203485 M. : 9979858530

Works & Display :

C/127, Vishala Industrial Estate(E),
Outside Kathawada G.I.D.C., S.P. Ring Road,
Odhav, Ahmedabad-342415. Ph : +91-79-29092177

E-mail : accumaxsurgicals@hotmail.com

Website : www.accumaxsurgicals.com

Products and specifications are Subject to change without notice