

# LABOR TABLE

## TWO SECTION LABOR TABLE

**TSLT SS 01** - Two Section Labor Table - Fully SS with Mattress & Padded legrest

**TSLT SS TP 01** - Two Section Labor Table - with Superimposed SS Top

**TSLT MS 01** - Two Section Labor Table - Fully MSPC



- Available in Fully S.S. Structure as well as M S P C with Superimposed S.S. Top
- New improved design with Superb finish and excellent workmanship
- Whole table top capable of being tilted by suitable screw mechanism
- Foot section capable of being folded underneath to suit Lithotomy position
- Back Section adjustable on rack
- Protective side rails on three sides
- Full Size : 72" x 30" x 30"
- Available with Mattress and Padded Leg Rest (Only in Fully S.S., Not in S.S. Top Model)

## PLAIN LABOR TABLE

PL LT 01



- M.S. Powder Coated Structure
- M S C P Lithotomy Rod with Straps with Perineal cut on leg Side
- Over all Size : 72" x 30" x 32"
- Optionally Available : Superimposed S.S. Top & Mattress

## ACCUMAX SURGICALS

Manufacturer of : HOSPITAL EQUIPMENTS & FURNITURE



**ACCUMAX**

### Admin & Works :

55, Amar Estate, Opp. Forge & Blower Comp,  
Nr. Memco, Naroda Road, Ahmedabad-25. INDIA.  
Ph : +91-79-22203485

### Works & Display :

C/127, Vishala Ind. Estate,  
Outside Kathawada G.I.D.C., S.P Ring Road,  
Odhav, Ahmedabad-342415. Ph : +91-7929092177

Mo.: +91 99798 58530 | E-mail : [accumaxsurgicalsh@hotmail.com](mailto:accumaxsurgicalsh@hotmail.com) | Website : [www.accumaxsurgicalsh.com](http://www.accumaxsurgicalsh.com)